
- EN:** Instructions for Use  
Ankura™ TAA Stent Graft System
- BG:** Инструкции за употреба за  
стент-графт система Ankura™ TAA
- CS:** Návod k použití  
systému Ankura™ TAA stent-graft
- DA:** Brugsanvisning  
Ankura™ TAA-stenttransplantatsystemer
- DE:** Gebrauchsanweisung  
Ankura™ TAA-Stentgraftsyste
- EL:** Οδηγίες χρήσης  
Συστήματα ενδομοσχευμάτων Ankura™ TAA
- ES:** Instrucciones de uso  
Sistemas de endoprótesis TAA Ankura™
- FR:** Notice d'utilisation  
Systèmes d'endoprothèse Ankura™ TAA
- HR:** Upute za uporabu  
Ankura™ TAA Stent Graft sustavi
- IT:** Istruzioni per l'uso  
Ankura™ Sistemi di endoprotesi TAA
- NL:** Gebruiksaanwijzing  
Ankura™ TAA stenttransplantaatsystemen
- PL:** Instrukcja obsługi  
Systemy stentgraftu TAA Ankura™
- PT:** Instruções de utilização  
Ankura™ Sistemas de Endoprótese TAA
- RO:** Instrucțiuni de utilizare  
Sisteme de greță stent TAA Ankura™
- SV:** Bruksanvisning  
Stentgraftsystemet Ankura™ TAA
- TR:** Kullanım Talimatları  
Ankura™ TAA Stent Greft Sistemleri

	<p>Catalogue number   Каталоген номер   Katalogové číslo   Katalognummer   Katalognummer   Αριθμός καταλόγου   Número de referencia   Référence   Kataloški broj   Numero di catalogo   Catalogusnummer   Numer katalogowy   Número do catálogo   Număr catalog   Artikelnummer   Katalog numarası</p>
	<p>Use-by date   Срок на годност   Datum spotřeby   Sidste anvendelsesdato   Verfallsdatum   Ημερομηνία λήξης   Fecha de caducidad   Date de péremption   Rok upotrebe   Data di scadenza   Houdbaarheidsdatum   Termin ważności   Prazo de validade   A se utiliza înainte de   Förbrukningsdatum   Son kullanma tarihi</p>
	<p>Date of manufacture   Дата на производство   Datum výroby   Produktionsdatum   Herstellungsdatum   Ημερομηνία κατασκευής   Fecha de fabricación   Date de fabrication   Datum proizvodnje   Data di fabbricazione   Productiedatum   Data produkcji   Data de fabrico   Data fabricației   Tillverkningsdatum   Üretim tarihi</p>
	<p>Serial number   Серийен номер   Sériové číslo   Seriennummer   Seriennummer   Αριθμός σειράς   Número de serie   Numéro de série   Serijski broj   Numero di serie   Seriennummer   Numer seryjny   Número de série   Număr de serie   Seriennummer   Seri numarası</p>
	<p>Lot Number   Номер на партида   Číslo šarže   Lot-nummer   Chargennummer   Αριθμός παρτίδας   Número de lote   Numéro de lot   Broj serije   Numero di lotto   Partijnummer   Numer partii   Número de lote   Număr lot   Batchnummer   Parti Numarası</p>
	<p>Consult Instructions for Use   Видже инструкциите за употреба   Přečtěte si návod k použití   Se brugsanvisningen   Gebrauchsanweisung beachten   Συμβουλευθείτε τις οδηγίες χρήσης   Consultar las instrucciones de uso   Lire la notice d'utilisation   Proučite Upute za uporabu   Consultare le istruzioni per l'uso   Raadpleeg de gebruiksaanwijzing   Patrz Instrukcja stosowania   Consultar as instruções de utilização   A se consulta instrucțiunile de utilizare   Se bruksanvisning   Kullanma Kılavuzuna Başvurun</p>

	<p>Double sterile barrier system   Система с двойна стерилна бариера   Systém dvojité sterilní bariéry   Dobbelt sterilt barrieresystem   System mit doppelter Sterilbarriere   Σύστημα διπλού αποστειρωμένου φραγμού   Sistema con doble barrera estéril   Système à double barrière stérile   Dvostruki sterilni sustav barijere   Sistema con doppia barriera sterile   Dubbel steriel barrièresysteem   System podwójnej bariery sterylnej   Sistema de barreira estéril dupla   Sistem de barieră sterilă dublă   Dubbla sterila barriärsystem   Çift steril bariyer sistemi</p>
	<p>Sterilized using ethylene oxide   Стерилизовано с этиленов оксид   Sterilizováno za použití ethylenoxidu   Steriliseret med ethylenoxid   Mit Ethylenoxid sterilisiert   Αποστειρωμένο με αιθυλενοξειδίο   Esterilizado con óxido de etileno   Stérilisé à l'oxyde d'éthylène   Sterilizirano etilen-oksidom   Sterilizzato con ossido di etilene   Gesteriliseerd met ethyleenoxide   Wyjałowiono tlenkiem etylenu   Esterilizado utilizando óxido de etileno   Sterilizat cu oxid de etilenă   Steriliserad med etylen   Etilen oksitle sterilize edilmıştır</p>
	<p>Non-pyrogenic   Непирогенно   Нерυρογενní   Ikke-pyrogen   Pyrogenfrei   Μη πυρετογόνο   Απιδόgeno   Non pyrogène   Nepirogeno   Απιδόgeno   Niet-pyrogeen   Produkt niepirogenny   Não pirogénico   Απιδόgen   Ikke-pyrogen   Pirojenik deǵildir</p>
	<p>Do not reuse   За еднократна употреба   Нерουžívejte opakovaně   Må ikke genbruges   Nicht wiederverwenden   Να μην επαναχρησιμοποιείται   No reutilizar   Ne pas réutiliser   Nije za višekratnu uporabu   Non riutilizzare   Niet opnieuw gebruiken   Nie używać ponownie   Não reutilizar   A nu se reutiliza   Återanvänd inte   Tekrar kullanmayın</p>

	<p>Do not use if package is damaged and consult instructions for use   Да не се използва, ако опаковката е повредена и се консултирайте с инструкциите за употреба   Pokud je balení poškozené, nepoužívejte produkt a přečtěte si návod k použití   Må ikke bruges, hvis emballagen er beskadiget. Se brugsanvisningen   Nicht verwenden, wenn die Verpackung beschädigt ist, und Gebrauchsanweisung beachten   Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί φθορά και συμβουλευθείτε τις οδηγίες χρήσης   No utilizar si el embalaje está dañado y consultar las instrucciones para su uso   Ne pas utiliser si l'emballage est endommagé et lire la notice d'utilisation   Nemojte koristiti ako je pakiranje oštećeno i pogledajte upute za uporabu   Non utilizzare se la confezione è danneggiata. Consultare le istruzioni per l'uso.   Niet gebruiken als de verpakking beschadigd is en raadpleeg de gebruiksaanwijzing   Nie używać, jeśli opakowanie jest uszkodzone i zapoznać się z instrukcją obsługi   Não utilizar se a embalagem estiver danificada e consultar as instruções de utilização   Nu utilizați dacă pachetul este avariât și consultați instrucțiunile de utilizare   Använd inte om förpackningen är skadad och se bruksanvisningen   Ambalaj hasar görmüşse kullanmayın ve kullanım talimatlarına başvurun</p>
	<p>Do not re-sterilize   Да не се стерилизира повторно   Neresterilizujte   Må ikke gensteriliseres   Nicht erneut sterilisieren   Να μην επαναποστειρώνεται   No volver a esterilizar   Ne pas restériliser   Ne višekratno sterilizirati   Non risterilizzare   Niet opnieuw steriliseren   Nie wyjalawiać   Não voltar a esterilizar   A nu se re-steriliza   Omsteriliserate inte   Tekrar sterilize etmeyin</p>
	<p>Store at room temperature   Да се съхранява на стайна температура   Skladujte při pokojové teplotě   Opbevares ved rumtemperatur   Bei Raumtemperatur aufbewahren   Φυλάσσεται σε θερμοκρασία δωματίου   Almacenar a temperatura ambiente   Stocker à température ambiante   Čuvati na sobnoj temperature   Conservare a temperatura ambiente   Bewaren op kamertemperatuur   Przechowywać w temperaturze pokojowej   Conservar à temperatura ambiente   A se păstra la temperatura camerei   Förvara i rumstemperatur   Oda sıcaklığında saklayın</p>

	<p>Keep dry   Да се държи сухо   Udržujte v suchu   Skal holdes tør   Trocken halten   Να διατηρείται στεγνό   Mantener seco   Conserver au sec   Čuvajte na suhom mjestu   Conservare in luogo asciutto   Droog bewaren   Przechowywać w suchym miejscu   Manter seco   A se pãstra ferit de umezealã   Förvaras torrt   Kuru ortamda saklayın</p>
	<p>Keep away from sunlight   Да се държи далеч от слънчева светлина   Uchovávejte mimo dosah slunečního záření   Må ikke udsættes for sollys   Von Sonnenlicht fernhalten   Να φυλάσσεται μακριά από το άμεσο ηλιακό φως   Mantener alejado de la luz solar   Conserver à l'abri de la lumière   Čuvajte dalje od sunčeve svjetlosti   Tenere lontano dalla luce del sole   Uit de buurt van zonlicht bewaren   Chronić przed promieniowaniem słonecznym   Manter afastado da luz solar   A se pãstra departe de lumina soarelui   Utsätt inte för solljus   Güneş ışığından koruyun</p>
	<p>Manufacturer   Производител   Výrobce   Producent   Hersteller   Κατασκευαστής   Fabricante   Fabricant   Proizvođač   Produttore   Fabrikant   Producent   Fabricante   Producător   Tillverkare   Üretici</p>
	<p>Authorized representative in the European Community   Оторизиран представител в Европейската общност   Oprávněný zástupce v Evropském společenství   Autoriseret repræsentant i det Europæiske Fællesskab   Autorisierte Vertretung in der Europäischen Gemeinschaft   Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Ένωση   Representante autorizado en la Unión Europea   Représentant autorisé dans la communauté européenne   Ovlašteni predstavnik u europskoj zajednici   Rappresentante autorizzato per la Comunità Europea   Bevoegde vertegenwoordiger in de Europese Gemeenschap   Autoryzowany przedstawiciel we Wspólnocie Europejskiej   Representante autorizado na Comunidade Europeia   Reprezentant autorizat în Comunitatea Europeană   Auktoriserad representant i Europeiska unionen   Avrupa Topluluğu'nda bulunan yetkili temsilci</p>

	<p>MR conditional   МР съвместим   Podmíněně bezpečný při magnetické rezonanci   MR-betinget   Bedingt MR-tauglich   Ασφαλές για μαγνητική τομογραφία υπό προϋποθέσεις (MR Conditional)   Condicional para RM   Compatible IRM sous conditions   Sigurno za uporabu u MR okruženju   A compatibilità RM condizionata   MR-condities   Warunkowo bezpieczny w środowisku MR   Condicionado por RM   Compatibilitate RM condiționată   MR-villkor   MR koşullu</p>
	<p>Made in China   Произведено в Китай   Vyrobeno v Číně   Fremstillet i Kina   Hergestellt in China   Κατασκευάζεται στην Κίνα   Fabricado en China   Fabriqué en Chine   Proizvedeno u Kini   Prodotto in Cina   Geproduceerd in China   Wyprodukowano w Chinach   Fabricado na China   Fabricat în China   Tillverkad i Kina   Çin mali</p>
	<p>Medical device   Медицинско изделие   Zdravotnický prostředek   Medicinsk udstyr   Medizinprodukt   Ιατροτεχνολογικό προϊόν   Producto sanitario   Dispositif médical   Medicinski uređaj   Dispositivo medico   Medisch instrument   Wyrób medyczny   Dispositivo médico   Dispozitiv medical   Medicinteknisk produkt   Tibbi cihaz</p>
	<p>Unique device identifier   Уникален идентификатор на изделието   Jedinečný identifikátor prostředku   Unik udstyrsidentifikation   Eindeutige Gerätekenung   Αποκλειστικό αναγνωριστικό τεχνολογικού προϊόντος   Identificador único de producto   Identifiant unique du dispositif   Jedinstveni identifikator uređaja   Identificazione unica del dispositivo   Unieke apparaat-ID   Unikalny identyfikator urządzenia   Identificador único do dispositivo   Identicator unic dispozitiv   Unik enhetsidentifikation   Benzersiz cihaz tanımlayıcısı</p>
	<p>Importer   Вносител   Dovozece   Importør   Importeur   Εισαγωγέας   Importador   Importateur   Uvoznik   Importatore   Importeur   Importer   Importador   Importator   Importör   İthalatçı</p>

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
Read all instructions carefully. Failure to follow the instructions, warnings and precautions may lead to serious complications or injury to the patients.

**1. Device Description**

- Ankura™ TAA Stent Graft System
  - ◆ The Ankura™ TAA Stent Graft System includes:
    - ◆ The Ankura™ TAA Stent Graft
    - ◆ The Ankura™ Delivery System
  - ◆ The Ankura™ TAA Stent Graft is pre-loaded into the Delivery System. The loaded delivery system is inserted endoluminally via the femoral or iliac artery and tracked through the patient's vasculature to deliver the stent graft to the target site.
- Ankura™ TAA Stent Graft


**Fig.1 Ankura™ TAA Stent Graft**

**Table 1 Stent Graft Materials**

Component	Material	Concentration in wt. %
Support Spring	Nitinol	49-67
Stainless steel sleeve	316LVM	7-17
Graft	PTFE	20-40
Marker	Platinum-Iridium	0.5-6
Suture	Polypropylene	<1


- ◆ The Ankura™ stent graft is a tube composed of ePTFE film supported by a metal mesh called a stent. The markers are fixed on the stent. The graft which is ultrathin artificial vascular material fixed to the stent without suture, can isolate aneurysm from blood, and it is stronger than the weakened artery and it allows the blood to pass through it without pushing on the bulge. The stent can provide stable and expansionary force necessary to open the graft.
- Ankura™ Delivery System

1	Tapered Tip	2	Markband	3	Proximal Anchor
4	Stent Graft	5	Inner Tube	6	Immobile Pole
7	Sheath	8	Front Grip	9	Slider Grip
10	Trigger	11	Screw Gear	12	Haemostatic Vale
13	Rear Grip	14	Proximal Releaser	15	Secure Buckle
16	Luer Connector				

**Fig.2 Ankura™ Delivery System**

- ◆ The Ankura™ stent graft is loaded inside the Ankura™ Delivery System. The Ankura™ Delivery System facilitates the placement of the stent graft via the arterial vasculature (e.g., femoral arteries). Using fluoroscopic guidance, the Ankura™ Delivery System is properly positioned within the patient's vasculature and the stent graft is deployed from the delivery system.
- Device model and specification
  - ◆ The model and specification of Ankura™ TAA Stent Graft System is indicated as follow and the details are shown in table 2.

**Table2 Summary of Ankura™ TAA Stent Graft System**

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
20	20	40,60,80,100,120	20F
22	22	40,60,80,100,120,	20F
24	20	60,80,100,120,140,160	20F
	24	40,60,80,100,120,140,160	
26	22	60,80,100,120,140,160	20F
	26	40,60,80,100,120,140,160	
28	22	120,140,160,180,200	20F
	24	60,80,100,120,140,160,180,200	
	28	40,60,80,100,120,140,160,180,200	
30	22	120,140,160,180,200	20F
	24	120,140,160,180,200	
	26	60,80,100,120,140,160,180,200	
	30	40,60,80,100,120,140,160,180,200	
32	24	120,140,160,180,200	20F
	26	120,140,160,180,200	
	28	120,140,160,180,200	
	32	40,60,80,100,120,140,160,180,200	

Proximal Diameter (D0/mm)	Distal Diameter (D1/mm)	Covered Length (L1/mm)	Delivery System
34	26	120,140,160,180,200	20F
	28	120,140,160,180,200	
	30	120,140,160,180,200	
	34	40,60,80,100,120,140,160,180,200	
36	28	120,140,160,180,200	20F
	30	120,140,160,180,200	
	32	120,140,160,180,200	
	36	40,60,80,100,120,140,160,180,200	
38	30	120,140,160,180,200	22F
	32	120,140,160,180,200	
	34	120,140,160,180,200	
	38	40,60,80,100,120,140,160,180,200	
40	32	120,140,160,180,200	22F
	34	120,140,160,180,200	
	36	120,140,160,180,200	
	40	40,60,80,100,120,140,160,180,200	
42	34	120,140,160,180,200	22F
	36	120,140,160,180,200	
	38	120,140,160,180,200	
	42	40,60,80,100,120,140,160,180,200	
44	36	120,140,160,180,200	22F
	38	120,140,160,180,200	
	40	120,140,160,180,200	
	44	40,60,80,100,120,140,160,180,200	
46	38	120,140,160,180,200	22F
	40	120,140,160,180,200	
	42	120,140,160,180,200	
	46	40,60,80,100,120,140,160,180,200	

Physicians should strictly adhere to indications and contraindications.

## 2. Indications for Use

- Ankura™ TAA Stent Graft System is indicated for endovascular repair of patients with thoracic aortic aneurysms or thoracic aortic dissection, requiring:
  - ◆ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
    - ◆ Aortic inner diameter in the range of 18-44mm.
    - ◆  $\geq 15$ mm non-aneurysmal aorta proximal and distal to the lesion.
    - ◆ Morphology suitable for endovascular repair.
  - ◆ Type B dissections in patients who have appropriate anatomy, including;
    - ◆ Adequate iliac/femoral access vessel that is compatible with the required delivery system.
    - ◆  $\geq 15$ mm landing zone proximal to the primary entry tear; proximal extent of the landing zone must not be dissected. Diameter at proximal extent of proximal landing zone in the range of 16-44mm.
    - ◆  $\geq 15$ mm landing zone distal to the primary entry

tear; distal extent of the landing zone must not be dissected. Diameter at distal extent of distal landing zone in the range of 16-44mm.

- ◆ Morphology suitable for endovascular repair.
- Additional considerations for patient selection include:
  - ◆ Patient's age and life expectancy;
  - ◆ Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity);
  - ◆ Toleration to the general, regional or local anesthesia;

The final decision is at the discretion of the physician and patient.

- Device Selection
  - ◆ The size of the Ankura™ Stent Graft should be appropriate to fit the patient's aortic anatomy. For aneurysms, appropriate oversizing (10%-20%) is built into the recommended sizes. For Type B dissections, appropriate oversizing (5%-20%) is built into the recommended sizes.

### 3. Contraindications

- The Ankura™ TAA Stent Graft is contraindicated in:
  - ◆ Patients with acute systemic infection.
  - ◆ Patients who have had other devices implanted in the cardiovascular cavity, which will interfere with the placement of this device.
  - ◆ Patients with mesenteric blood flow mainly supplied by the inferior mesenteric artery;
  - ◆ Patients who have allergic reaction to the device;
  - ◆ Patients who are not suitable for endovascular repair in vascular morphology;
  - ◆ Patients who can not tolerate contrast agents due to renal insufficiency;
  - ◆ Patients who are allergic to contrast agents;
  - ◆ Aneurysms neck with thrombus.
  - ◆ Non-aneurysmal aortic proximal neck length <1.5cm.
  - ◆ Non-aneurysmal aortic distal anchorage zone <1.5cm
  - ◆ For aneurysms, Non-aneurysmal aortic diameter <18mm or >42mm. For Type B dissections, Non-aneurysmal aortic diameter <16mm or >44mm.

### 4. Pre-operative Preparation

Read the instructions carefully before using the device. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help to guide the physician and do not replace the physician judgment.

- Materials required
  - ◆ 0.035"/260cm Super Stiff Guidewire. For example, Lifetech's Super Stiff Guidewire, Amplatz Super Stiff Guidewire or comparable.
  - ◆ Fluoroscope with digital angiography capabilities;
  - ◆ Angiographic catheter, and marker catheter;
  - ◆ Contrast agent;
  - ◆ High pressure Syringe;
  - ◆ Heparinized saline solution;
  - ◆ Any other devices required by the physician for the operation (including devices used for femoral artery dissection).
- A multidisciplinary team that has combined procedural experience with:
  - ◆ Local anesthesia or general anesthesia;
  - ◆ Femoral artery exposure, arteriotomy and repair;
  - ◆ Percutaneous access and closure techniques;
  - ◆ Fluoroscopic and angiographic image interpretation;
  - ◆ Endovascular stent placement;
  - ◆ Appropriate use of radiographic contrast agent;
  - ◆ Techniques to minimize radiation exposure;
  - ◆ Expertise in necessary patient follow-up modalities;
  - ◆ Experience in traditional vascular surgery for any possibility of failure in EVAR.
- Products preparation
  - ◆ The Stent Graft diameter should be selected preferred to vessel inner diameter not outer diameter. The uncorrectable size of the device may result in endoleak or stent-graft migration. Generally, the diameter of Stent Graft selected should be 1.1 to 1.2 times for aneurysm,

1.05 to 1.2 times for type B dissections times larger than the landing zone vessel inner diameter. Recommended reference sizes for aneurysms and aortic dissections are shown in table 3.

**Table 3 Recommended reference sizes for aneurysms and aortic dissections**


Recommended reference size of aneurysm		Recommended reference size for aortic dissection	
Vessel diameter in anchoring area (mm)	Diameter of stent graft (mm)	Vessel diameter in anchoring area (mm)	Diameter of stent graft (mm)
17-18	20	17-19	20
18-20	22	19-21	22
20-22	24	21-23	24
22-24	26	23-25	26
24-25	28	24-27	28
25-27	30	26-29	30
27-29	32	28-30	32
28-31	34	30-32	34
30-33	36	31-34	36
32-35	38	33-36	38
33-36	40	35-38	40
35-38	42	37-40	42
37-40	44	38-42	44
38-42	46	40-44	46

Note: the attached table is only the recommended reference size, and the final selection of stent graft specifications is subject to the doctor's judgment based on the actual situation of the patient.

- ◆ Inspect the device and the package to verify that no damage has occurred as a result of shipping. Do not use the device if any damage of the device or the sterilization barrier has been observed. If damage has occurred, return it to Lifetech. Prior to use the device, re-review the specification conforming to the patient required.
- Patient preparation
  - ◆ Measure proximal aortic neck length and vessel inner diameter. Aortic neck diameters should be measured from axial CTA films. For aneurysms, two diameter measurements are required for both the proximal and distal necks. All measurements per neck must be within one intended aortic inner diameter range. For Type B dissections, two diameter measurements are required for both the proximal extent of the proximal neck and the distal extent of the distal which must be in non-dissected tissue.
  - ◆ Measure the angle of the proximal aortic neck relative to the long axis of aneurysm.
  - ◆ Evaluate the quality of aortic neck.
  - ◆ Assess the morphology of the vessel compatible with endovascular aortic repair.
  - ◆ Undergo any necessary common examination, and monitoring of vital signs.
  - ◆ Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
  - ◆ Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral artery bifurcations.
  - ◆ Adequate iliac/femoral artery access and expose one common femoral artery using standard surgical technique.
  - ◆ Systemic anticoagulation should be used during the deployment procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.


## 5. Operation Procedure

- Preparation before implant
  - ◆ Attach the syringe with heparinized saline to the luer hub and hemostatic valve and flush separately to evacuate air from the inner tube and the outer sheath.(Fig.3)


**Fig 3 Flush the inner tube and the outer sheath**

- ◆ Monitor important vital signs.
- ◆ Administer heparin to adjust ACT (activated clotting time) to be 200 seconds. Inspect flush solution; sheath and guidewire should be flushed after each exchange.
- ◆ Use gauze pads with saline solution to wipe the outer sheath and activate the hydrophilic coating. (Fig. 4)


**Fig. 4 Wipe the sheath with sterile heparinized normal saline**

- ◆ Expose one side femoral arteries by standard surgery.
- ◆ Place a 0.035"/260cm super stiff guidewire in the femoral artery and advance it until to the ascending artery. (Fig. 5)

**Fig. 5 The super stiff guide wire enters the predetermined position**

- Deploy the stent graft
  - ◆ Advance the delivery system to the aneurysms or dissection sites along the guidewire, and always maintain the hemostatic valve toward the left side of the patient (Fig. 6)


**Fig. 6 Introduction of the delivery system**

**Notes:**

- ◆ The connecting bar of the TAA Stent Graft must be in the outer curved side of the aortic arch. And the connecting bar has the same side with the hemostatic valve of the delivery system when the TAA Stent Graft is preloaded into the delivery system.
- ◆ Due to vascular distortion, the hemostatic valve of the delivery system may not be completely on the same side as the connecting bar of the Stent Graft. In case of this situation, the positioning of 8-shaped marker on the stent graft shall be used to ensure that the connecting bar of the TAA Stent Graft is locating in the outer curved


side of the aortic arch.

- ◆ Ensure the Stent Graft positioning at the aneurysms or dissection sites. Ensure the connecting bar of the TAA Stent Graft in the outer curved side of aortic arch. Confirmation method: the 8-shaped marker at the proximal end of the stent graft is located in the outer curved side of the aortic arch.
- ◆ Repeat the angiography to verify the delivery system locating at the aneurysms or dissection sites.
- ◆ Decrease the systolic blood pressure (SBP) to 10.5 to 12 kpa.
- ◆ Slowly deploy the stent graft. Hold the front grip of the delivery system stationary, then rotate the slider grip counterclockwise to slowly deploy the proximal part of the stent graft until the covered stent adhere to the vascular wall. (Fig. 7)


**Fig. 7 Slowly deploy the stent graft**

- ◆ Rapidly deploy the stent graft. Hold the front grip of the delivery system stationary, then unlock the trigger and pull the slider grip along the screw gear to fully deploy the stent graft. (Fig. 8)


**Fig. 8 Rapidly deploy the stent graft**

- Proximal end release
  - ◆ Keep the delivery system stationary. Remove the Secure Buckle from the grip of the delivery system. (Fig. 9)

**Fig. 9 Remove the Secure Buckle**

- ◆ Pull the Proximal Releaser to release the proximal end of the stent graft. (Fig. 10)


**Fig. 10 Pull the Proximal Releaser**

- ◆ Verify the position and the state of the completely deployed stent graft under fluoroscopy.

- Remove the delivery system
  - ◆ Hold the slider grip of the delivery system stationary, unlock the trigger and pull back the front grip to make the tip back to the distal end of the out sheath. (Fig.11)

**Fig. 11 Pull back the tip**

**Note: Pull back the front grip should be gently and under fluoroscopy guidance to avoid stent graft migration.**

- ◆ Remove the delivery system and guidewire completely and close the puncture site with appropriate techniques.
- Handle Disassembly Technique for Partial Stent Graft Deployment

In the unlikely event of delivery system failure and concomitant partial stent graft deployment, a “handle disassembly” technique may permit the successful deployment of the stent graft.

Step 1 - Pull back the trigger and fully retract the slider. Note: Since the graft cover is severed, the slider can be retracted without further deploying the stent graft.

Step 2 - Stabilize the delivery system.

Step 3 - Insert the tips of a pair of hemostats into gap at the end handle on the slider grip. Then, slide the hemostatic forward to split the slider grip into left and right halves.

Step 4 - Remove the left and right halves of the slider grip from the screw gear to expose the sheath joint.

Step 5 – Hold the front handle with one hand and hold the sheath joint with the other hand and pull back until the stent graft is fully deployed.

Step 6 – Remove the delivery system by gripping the screw gear and withdrawing from the patient.

## 6. Potential Adverse Events

Adverse events associated with use of the Ankura™ TAA Stent Graft System include, but are not limited to the following:

- Allergic reaction
- Aneurysms enlargement
- Aneurysms/vessel rupture
- Anesthetic complications
- Arterial or venous thrombosis (including thrombosis and embolization)
- Arterial damage
- Arterial perforation
- Arterial stenosis
- Arrhythmia
- Arteriovenous fistula or pseudoaneurysm
- Bleeding/Hematoma
- Bowel complications
- Branch vessel occlusion
- Cardiac failure/Infarction
- Coagulopathy
- Conversion to open surgery
- Death
- Edema
- Embolism
- Hemorrhage/bleeding
- Hypotension / hypertension
- Infection
- Intercostal pain
- New dissecting aneurysm
- Neurological damage, local or systemic (e.g., stroke, paraplegia, paraparesis)
- Pulmonary complications
- Paraparesis
- Paraplegia
- Paralysis

- Paresthesia
- Pericardial effusion
- Renal function insufficiency/kidney failure
- Retrograde Type A dissection
- Change in mental status
- Reverse tearing or rupture of dissection
- Shock
- Tissue necrosis
- Wound complications and subsequent attendant problems
- Device complications due to:
  - ◆ Insertion and removal difficulties
  - ◆ Release difficulties
  - ◆ Failure to deploy
  - ◆ Inaccurate placement
  - ◆ Endoleak
  - ◆ Stent migration
  - ◆ Stent fracture
  - ◆ Rupture of graft material
  - ◆ Erosion
  - ◆ Stent thrombosis or embolization
  - ◆ Stent graft expansion
  - ◆ Stent graft twisting or kinking

**Note:**

For any serious incident, the user or patient should report it to Lifetech at [quality@lifetechmed.com](mailto:quality@lifetechmed.com) or the local sale representative of Lifetech and your local Competent Authority (CA).

**7. Summary of Clinical Data**

Three clinical studies of Ankura™ TAA stent graft system have been finished in China, including a mid-long-term (one-to-nine-year) follow-up study. These studies aimed to evaluate the clinical safety and efficacy of the Ankura™ Stent Graft Systems for endovascular treatment of thoracic aortic aneurysm and Stanford B-type dissection, and they met the predefined primary endpoints of treatment success with few adverse events.

**Study 1**

107 patients implanted Ankura™ stent graft systems (TAA-62 cases, AAA-30 cases, AUJ-15 cases) were enrolled in the study in China. The stent grafts were successfully implanted in all 107 patients (100%). 91 of 106 patients were followed up after the procedure, and the mean follow-up time was 218 (range 30 – 600) days. Treatment success rate of aortic aneurysm at 6-month post-implant reached 93.4% (85/91). 85 of 91 patients evaluated for safety events during the follow up were free of a safety event (93.4%).

Follow-up data showed that:

One out of 107 patients (0.9%) died in-hospital (0.9%), but this death was unrelated to the device or procedure.

One out of 107 patients (0.9%) had artery stenosis.

One out of 107 patients (0.9%) experienced paraplegia.

One out of 107 patients (0.9%) had renal failure.

Three out of 107 patients (2.8%) experienced early endoleaks.

No patient experienced stent displacement, stent fracture, or aortic rupture.

**Study 2**

This was a randomized control trial. A total of 144 patients with aortic aneurysm or dissection underwent endovascular aneurysm treatment with implantation of stent grafts. 71 (49.3%) patients were implanted with Ankura™ TAA stent grafts. The immediate technical success rate was 100%. Of the 65 patients implanted Ankura™ TAA stent grafts evaluated for safety events at one year, 58 were free of a safety event (90%).

One-year patient follow-up has detected the following adverse events: respiratory failure and heart failure (N=1), paraplegia (N=1), death (N=3), pulmonary infection (N=1), abdominal distension (N=1). One out of 71 patients implanted the Ankura™ TAA stent experienced a device-related adverse event. CTA follow-up showed no displacement, stent thrombosis, stent discounts, stent rupture, and the rate of internal leakage.

**Study 3**

This study retrospectively collected mid-long-term real-world clinical data from multiple hospitals in China. 305 patients were enrolled and they completed the mid- and long-term follow-up.

The immediate technical success rate was 98.69% (301/305). Thoracic aortic dissection was successfully treated in 110 cases 2 years after operation, with a mid-long-term success rate of 94.02% (110/116).

179 of 305 patients have one-year or more follow-up data, and the data showed that:

Two out of 305 patients (0.66%) experienced aortic dissection-related death within 30 days after operation.

Two out of 305 patients (0.66%) experienced cerebral infarction within 30 days after operation.

One out of 305 patients (0.33%) experienced paraplegia within 30 days after operation.

One out of 305 patients (0.33%) experienced thoracic aortic rupture.

There were 19 deaths (6.23%) in this study, 17 of which were not device-related. For the other two deaths, there was insufficient evidence to judge their relationship with the device.

No patient experience stent migration, stent fracture, or type I/III endoleaks in the mid-long-term (one year and more).

## **8. Storage and Transport Requirements**

- The Ankura™ TAA Stent Graft is pre-loaded in delivery system, and stored in sterile package.
- This device should be stored in a clean, well ventilated room without corrosive gas and stored at room temperature.
- During transportation, it is strictly prohibited to squeeze, get wet and expose to the sun, and handle with care.
- No exposure to organic solvent.

## **9. Warnings**

- General
  - ◆ Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
  - ◆ The Ankura™ TAA Stent Graft can not be used after resterilization.
  - ◆ If any sign of damage is observed prior to use, do not use the product and return it to Lifetech.
  - ◆ This device should only be used by physicians and teams trained in vascular interventional techniques and in the use of device.
  - ◆ The long term performance of endovascular grafts has not yet been confirmed, all patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft (e.g., endoleak, migration, distortion, aneurysms enlarging dissection tearing etc.).
  - ◆ Clinical studies show that uncorrectable patient selection may result in poor curative effect. In addition, if unsuccessful endovascular aortic repair require open surgery may result in an increased risk of complications relative to the two procedures.
  - ◆ Patient's condition needs to be evaluated before operation on the following aspects, cardiac failure/infarction, pulmonary complications, edema, and renal function insufficiency.
  - ◆ After the device placement, patients should be regularly monitored for the circumstances of the device (including endoleak, migration, aneurysms enlarging dissection tearing etc.). At a minimum, annual X-ray, contrast or non-contrast CT imaging is required.
  - ◆ Patients with specific clinical findings, e.g, enlarging aneurysms or tearing dissection or unexpected diseases in the position of the endovascular graft should undergo secondary interventions or open surgical treatment. An increase in aneurysms or tearing dissection size and/or persistent endoleak may lead to aneurysms or tearing dissection rupture.
  - ◆ Patients with infection after the device placement should administer oral antibiotics.
  - ◆ Always have a vascular surgery team available during stent-graft implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.
  - ◆ Do not re-use or re-sterilize the Ankura™ TAA Stent Graft System.

- Patients selection, treatment and follow-up
  - ◆ The safety and effectiveness of the Ankura™ TAA Stent Graft has not yet been evaluated in following populations:
    - ◆ Traumatic aortic injury
    - ◆ Uncorrectable coagulopathy
    - ◆ Hereditary connective tissue disease (e.g., Marfan's or Ehlers-Danlos Syndromes)
    - ◆ Patients with active systemic infections
    - ◆ Pregnant or nursing females
    - ◆ Morbid obesity
    - ◆ Less than 18 years of age
    - ◆ Patient's life expectancy less than one year
  - ◆ Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and the profile of the delivery system. Too tortuous vasculature may increase the failure in stent graft deployment or proximal end release.
  - ◆ The Ankura™ TAA Stent Graft System is not recommended in patients who weight greater than 350 Pounds (150Kg) or can not undergo accurate fluoroscopy examination due to obesity.
  - ◆ Patients with a systemic infection may be at increased risk of endovascular stent-graft infection.
  - ◆ Patients with uncorrectable coagulopathy may have an increased risk of Type II endoleak or bleeding complications.
  - ◆ The blood pressure of patients shall be measured every day after surgery, and clinical record shall be taken. High blood pressure may lead to stent expansion, vascular rupture and other complications. Blood pressure should be monitored in time for the medical staff to know the postoperative status of the patient. At the same time, the doctor needs to remind the patient to control blood pressure in the normal range.
  - ◆ The product should not be used in patients who are allergic to the materials listed out in Table 1.
- Implant Procedure
  - ◆ Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
  - ◆ Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
  - ◆ Maintain the guide wire position during delivery system insertion.
  - ◆ Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the stent graft.
  - ◆ Always use fluoroscopy for guidance, delivery and observation of the device within the vasculature.
  - ◆ The use of the Ankura™ TAA Stent Graft System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast agent used during the procedure.
  - ◆ To avoid any twist in the endovascular stent graft, during any rotation of the delivery system, be careful to rotate all the components of the system together.
  - ◆ Inaccuracy placement of the stent graft within the vessel may result in increased risk of endoleak, migration, or embolization.
  - ◆ Inadequate fixation of the stent graft may result in increased risk of migration of the stent graft, and incorrect deployment of the stent graft requires surgical intervention.
  - ◆ Do not continue advancing if resistance is felt during advancement of the guide wire or delivery system. Stop and assess the cause of resistance.
  - ◆ Unless medically indicated, do not deploy the stent graft in a location that will occlude arteries necessary to

supply blood flow to organs or extremities.

- **MRI Information**

Nonclinical testing has demonstrated that the Ankura™ Stent Graft is MR Conditional. A patient implanted with this device can be scanned safely under the following conditions:

- ◆ Static magnetic field of 3.0 Tesla or less
- ◆ Spatial gradient field of  $\leq 720$  Gauss/cm
- ◆ Maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning

Based on non-clinical testing, the device was determined to produce a temperature rise of less than 2°C at a maximum whole-body-averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of MR scanning. The maximum whole-body-averaged specific absorption rate (SAR) was derived by calculation and verified by calorimetry.

**Note:** MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant. The image artifact extends approximately 20 mm outside the device lumen when scanned in non-clinical testing using the spin echo sequence and gradient echo sequence, respectively in a 3.0 T Scanner with a whole body coil.

#### **10. Patient Counseling Information**

The biggest benefits of treatment of the patient's diseased or injured aorta with Ankura™ stent graft are decreased chance of rupture and restoration of normal blood flow. If left untreated, aortic lesions can expand and rupture, resulting in bleeding inside the body, which is life-threatening.

The physician should consider, but not be limited to, the following when explaining this endovascular device and related procedures to the patient:

- The differences between endovascular repair and open surgical repair, include the risks of open surgical repair and endovascular repair.
- Pros and cons of endovascular repair and open surgery repair.
- Endovascular repair has the potential advantage of minimal trauma.
- Endovascular repair or open surgical repair of the lesion may be required in the future.
- Details of possible complications after implantation of the device.
- Regular follow-up after implantation of the device.

#### **11. Follow Up**

- Periods of follow up
  - ◇ Clinical follow up should be performed at 24 hours, 1, 3, 6, 12 and 24 months after the procedure. The follow up period can be adjusted at any time according to the circumstances of each individual patient.
- Methods used during follow up
  - ◇ Haematological examination, X-ray, B-US, CTA, MRI, DSA etc.
- Follow up items:
  - ◆ Size of aneurysms or dissection;
  - ◆ Embolization;
  - ◆ Changing of aneurysms pulsatility or dissection;
  - ◆ Migration;
  - ◆ Endoleak;
  - ◆ Stent Graft distortion.
- ◇ In addition, patients with following specific findings during following up, endovascular intervention or surgical treatment should be considered:
  - ◆ Aneurysms enlarge 5mm and above;
  - ◆ Changing of aneurysms or dissection pulsatility with or without aneurysms enlarging or endoleak or dissection tearing;
  - ◆ Continual endoleak with or without aneurysms enlarging or dissection tearing;
  - ◆ Stent Graft migration results in endoleak.

**12. Shelf Life**

- The Ankura™ TAA Stent Graft System is sterilized with ethylene oxide. The shelf life is indicated on the label. Do not use the expired product.

**13. Package and Label**

- The Ankura™ TAA Stent Graft System is supplied sterile. Each Ankura™ TAA Stent Graft is individually contained within a Delivery System.
- The Ankura™ TAA Stent Graft with Delivery System is immobilized and protected in a PETG tray, then sealed with two Tyvek1073B pouches, on which a primary label is stuck. The product is sterilized and put in a box with IFU, implant card, customer feedback form and other accompanying document. A label is stuck on the box.

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Lifetech Scientific (Shenzhen) Co., Ltd